

Panreac

DERQUIM

Laboratory glassware cleaning agents



DERQUIM products have been specially made for cleaning the laboratory equipment. They can be used in the manual cleaning (series LM) or for washing using automatic machines (series LA). The various formulations are specifically adapted to any kind of laboratory residue: chemical, biological or clinical.

Two products with disinfecting and bactericidal properties are included in our DERQUIM program. DERQUIM DSF 01 is used for hand disinfecting and DERQUIM DSF 11 is used for surfaces and instruments disinfecting.

You can find detailed specifications of these products in our General Catalogue or in www.panreac.com

The following table shows our DERQUIM program.



Code	Description	Package					
		1 l	2 kg	3 l	5 l	10 kg/l	25 kg/l
For automatic or manual cleaning							
503574	DERQUIM + Universal Detergent, LIQUID						
For manual cleaning							
502600	DERQUIM LM 01 Alkaline LIQUID						
502601	DERQUIM LM 02 Neutral, phosphates free LIQUID						
502602	DERQUIM LM 03 Phosphates free LIQUID						
For automatic cleaning							
502603	DERQUIM LA 11 Slightly alkaline SOLID						
502604	DERQUIM LA 12 Alkaline SOLID						
502605	DERQUIM LA 13 Alkaline with detergents SOLID						
502606	DERQUIM LA 14 Slightly alkaline LIQUID						
502607	DERQUIM LA 15 Alkaline LIQUID						
502608	DERQUIM LA 21 Acid, with phosphoric acid LIQUID						
502609	DERQUIM LA 22 Acid, with citric acid LIQUID						
Automatic cleaning aids							
502610	DERQUIM LA 31 Antifoaming						
502611	DERQUIM LA 32 Rinsing						
503468	DERQUIM SALT (Sodium Chloride lumps)						
Disinfectants							
504993	DERQUIM DSF 01 Antiseptic for hands LIQUID						
504994	DERQUIM DSF 11 Antiseptic for surfaces and implements LIQUID						
Chromic Mixture							
502612	DERQUIM MC Chromic mixture						
506071	DERQUIM OXY (Substitute of Chromic mixture) SOLID						12 x 10 g

LM : Manual cleaning

LA : Automatic cleaning

DSF : Disinfectant

Packaging symbols:

Polyethylene bottle

Polyethylene canister

Glass bottle coated with PE

Polypropylene bucket with handle

Polyethylene container with removable tap

Box